

# Body Bag

## Standard



### SPECIFICATION

- 36" W x 90" L
- 7mm High Strength Virgin Vinyl
- Envelope Style Zipper for Easy Access
- #5 Nylon Zipper
- Includes: 3 Identification Tags

**Includes: 3 Identification Tags**



Item #	Description	UOM
MS-BOD100	Standard Body Bag	10/CS

# Body Bag

## Pediatric

**7mm Virgin Vinyl**  
High Strength



**Envelope Style Zipper**  
Easy Access

**Bound Edges**  
Extra Durability

**Includes: 3 Identification Tags**



### SPECIFICATION

- 33" W x 49" L
- 7mm High Strength Virgin Vinyl
- Envelope Style Zipper for Easy Access
- #5 Nylon Zipper
- Includes: 3 Identification Tags

Item #	Description	UOM
MS-BOD300	Pediatric Body Bag	10/CS