


Endotracheal Tubes

Uncuffed and Cuffed

MedSource Endotracheal Tubes are engineered with clear, flexible PVC, and a smooth Murphy eye as an alternative vent port



SPECIFICATION

- Standard 15mm connector
- Oral or nasal use
- Short or long-term use
- Medical grade PVC
- Sterilized for single-use
- Not made with natural rubber latex



www.medsourcelabs.com

Item #	ET Tube Uncuffed Size	UOM	Item #	ET Tube Cuffed Size
MS-23120	2.0mm	10 EA/BX 10 BX/CS 100 EA/CS	MS-23240	4.0mm
MS-23125	2.5mm		MS-23245	4.5mm
MS-23130	3.0mm		MS-23250	5.0mm
MS-23135	3.5mm		MS-23255	5.5mm
MS-23140	4.0mm		MS-23260	6.0mm
MS-23145	4.5mm		MS-23265	6.5mm
MS-23150	5.0mm	MS-23270	7.0mm	
MS-23155	5.5mm	MS-23275	7.5mm	
MS-23160	6.0mm	MS-23280	8.0mm	
		MS-23285	8.5mm	
		MS-23290	9.0mm	
		MS-23295	9.5mm	
		MS-23210	10.0mm	