



Barrier Sheet



5/8" hole in each corner for easy attachment to Blanket System



Highly absorbant interior holds fluids in place and prevents leaking

Fluid Resistant* Exterior keeps fluids safely locked inside

*(Blanket System not included)

The MedSource Barrier Sheet is a lightweight, fluid resistant* sheet that is manufactured to fit securely inside the Blanket System. The Barrier Sheet features an absorbant side made with hydrophillic PP for maximum absorbancy and a reverse side that is coated with fluid resistant* PE film. Each sheet is manufactured with 4 holes (5/8" diameter) for attaching to the blanket.

**Fluid could pass through the material with pressure or time*



www.medsourcelabs.com

Item #	Description	UOM
MS-42700	MedSource Barrier Sheet, 60" W x 83"H	50 EA/CS