


Total Lifter Patient Transport Pad



Provides Warmth
Reduces trauma impact

Absorbant Material
Up to 5.5 fl oz
(163 ml per square foot)



16 Non-Tear Handles
Supports up to 1000 lbs

SPECIFICATION

- 42" x 72"
- Material: Front: Polyester, Middle: Breathable Film, Bottom: Air-laid Paper
- 16 Non-Tear Handles
- Supports Up to 1000 lbs



www.medsourcelabs.com

Item #	Description	UOM
MS-95200	Total Lifter Patient Transport Pad	10 EA/CS