



SPECIFICATION

- Various sizes
- Malleable metal core coated in plastic
- Sterilized with ethylene oxide

Item #	Description	UOM
MS-ETS100	ET Stylet, 6FR	20 EA/BX, 10 BX/CS
MS-ETS200	ET Stylet, 10FR	20 EA/BX, 10 BX/CS
MS-ETS300	ET Stylet, 12FR	20 EA/BX, 10 BX/CS
MS-ETS400	ET Stylet, 14FR	20 EA/BX, 10 BX/CS