


# Penlight with Pupil Gauge



**Pupil Gauge**  
Provides an easy reference for the user



**Pocket Clip**  
Activates a concentric light beam

## SPECIFICATION

- 5" Long with 1/2" Diameter
- White with Blue Pupil Gauge
- Press the Pocket Clip for a Strong Concentric Light Beam



[www.medsourcelabs.com](http://www.medsourcelabs.com)

Item #	Description	UOM
MS-PEL100	Penlight with Eye Pupil Gauge	6 EA/PK, 50 PK/CS, 300 EA/CS