



One-Way Valve
Fits snugly into the re-usable body

Optional O2 Port
Allows for delivery of supplemental oxygen

Collapsible
Convenient to carry

SPECIFICATION

- Can Be Used on Infants, Children, and Adults
- Not Made With Natural Rubber Latex
- Acts As a Barrier Separating Rescuer and Patient
- Cuff is Pre-Inflated

Item #	Description	UOM
MS-PM103W	CPR Mask, White Case, With O2 Port	100 EA/CS
MS-PM104	CPR Mask, White Case, No O2 Port	100 EA/CS
MS-PMV	One Way Valve Filter Assembly	300 EA/CS