


VeroTrue®

Conventional IV Catheter



Available in PTFE or PUR Catheter Material

SPECIFICATIONS

- PTFE or PUR catheter material
- Ground stainless-steel needle
- Straight hub
- Disposable
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- ISO compliant color coding
- Radiopaque



www.medsourcelabs.com

VeroTrue® PTFE	Gauge and length	VeroTrue® PUR	Gauge and length	UOM
MS-85112	12G x 1.75"	MS-852142	14G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85113	13G x 1.75"	MS-852143	14G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85114	14G x 1.75"	MS-852162	16G x 1.16"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85116	16G x 1.75"	MS-85216	16G x 1.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85117	17G x 1.75"	MS-852164	16G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85118	18G x 1.75"	MS-85218	18G x 1.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85120	20G x 1"	MS-852184	18G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-851204	20G x 1.25"	MS-85220	20G x 1"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85122	22G x 1"	MS-852204	20G x 1.25"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85124	24G x 0.75"	MS-852202	20G x 1.88"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-85126	26G x 0.75"	MS-85222	22G x 1"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
		MS-852243	24G x 0.50"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
		MS-852242	24G x 0.56"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
		MS-85224	24G x 0.75"	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS