


# PurSafe

## Safety IV Cannula with Catheter and Without Injection Valve and Without Small wings\*



### SPECIFICATIONS

\*Not for sale in the U.S.

- PUR or PTFE catheter material
- Disposable
- Passive style safety mechanism
- Ground stainless-steel needle
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- Sterilized with ethylene oxide
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant
- Radiopaque

All catheters packaged 50/Box and 200 EA/CS

[www.medsourcelabs.com](http://www.medsourcelabs.com)

PUR MATERIAL			
Item #	Size	Item #	Size
MS-89114	14G x 1.25"	MS-891204	20G x 2.5"
MS-891142	14G x 1.77"	MS-891203	20G x 1.88"
MS-89116	16G x 1.25"	MS-891205	20G x 3"
MS-891162	16G x 1.77"	MS-89122	22G x 1"
MS-89118	18G x 1.77"	MS-891224	22G x 2.5"
MS-891182	18G x 1.25"	MS-891223	22G x 1.88"
MS-891184	18G x 2.5"	MS-891225	22G x 3"
MS-891183	18G x 3"	MS-89124	24G x 0.75"
MS-89120	20G x 1.25"	MS-891242	24G x 1.88"
MS-891202	20G x 1"	MS-891244	24G x 2.5"

PTFE MATERIAL			
Item #	Size	Item #	Size
MS-89314	14G x 1.25"	MS-893204	20G x 2.5"
MS-89314	14G x 1.77"	MS-893203	20G x 1.88"
MS-89316	16G x 1.25"	MS-893205	20G x 3"
MS-893162	16G x 1.77"	MS-89322	22G x 1"
MS-89318	18G x 1.77"	MS-893224	22G x 2.5"
MS-893182	18G x 1.25"	MS-893223	22G x 1.88"
MS-893184	18G x 2.5"	MS-893225	22G x 3"
MS-893183	18G x 3"	MS-89324	24G x 0.75"
MS-89320	20G x 1.25"	MS-893242	24G x 1.88"
MS-893202	20G x 1"	MS-893244	24G x 2.5"