


# TrueSafe Comfort<sup>®</sup> Neonatal



## SPECIFICATIONS

- Push button style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO compliant color-coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant



www.medsourcelabs.com

Item #	Gauge and length	Catheter ID (mm)	Catheter OD (mm)	Gravity flow rate (mL/min)	UOM
MS-84326	26G x 0.56 in	0.4 mm	0.6 mm	11 mL/min	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843262	26G x 0.75 in	0.4 mm	0.6 mm	11 mL/min	50 EA/BX, 4 BX/CS, 200 EA/CS