


Non-Contact Infrared Thermometer



Set Button
Choose between audible or silent mode

Color Digital Screen
Easy to understand

Memo Button
View up to 50 previous readings

Mode Button
Two temperature modes are available: Body and Surface

On/Scan Button
Easy one-handed operation

SPECIFICATION

- L 5.43" x W 3.74" x H 1.47"
- Weight: 3 Ounces (Without Batteries)
- Measurement Time: ≤ 2 Seconds
- Measurement Distance: 0.4" - 2.0"
- Body Temperature Range: 95.0° F - 107.6° F
- Ambient Temperature: 59° F - 95° F
- Accuracy: ± 0.4° F
- Automatic Power Off: ≤ 18 Seconds
- Requires 2 AAA Batteries

This is an FDA 510k cleared device.



www.medsourcelabs.com

Item #	Description	UOM
MS-131002	Non-Contact Infrared Thermometer, IR 300	40 EA/CS